visomat®

handy express

Die
mobile
Vorsorge-
Messung

Gebrauchsanweisung
Instructions for use
Mode d'emploi
Istruzioni per l'uso

DE-2 - DE-33
EN-34 - EN-63
FR-64 - FR-95
IT-96 - IT-127

**UEBE**
Germany est.1890

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Table of contents

| | | |
|----------|---|----|
| A | Intended use | 36 |
| B | Safety instructions | |
| | 1. Important patient information | 37 |
| | 2. Important technical information | 38 |
| C | Operating the device | |
| | 1. Device description | 40 |
| | 2. Display | 41 |
| | 3. Important usage information | 42 |
| | 4. Initial operation of the device | 43 |
| | 5. Inserting/replacing batteries, battery icon | 43 |
| | 6. Activating and setting the date/time | 43 |
| | 7. Attaching the cuff | 44 |
| | 8. Body position during the measurement | 45 |
| | 9. Selecting the user | 46 |
| | 10. Measuring blood pressure | 47 |
| | 11. Classification of the measured values (WHO) | 48 |
| | 12. Pulse pressure control | 49 |
| | 13. Cardiac arrhythmia/atrial fibrillation | 49 |
| | 14. Motion control | 50 |
| | 15. Using the memory | 51 |
| D | What you should know about blood pressure | |
| | 1. The systolic and diastolic blood pressure value | 53 |
| | 2. Why you measure different values | 53 |
| | 3. Why should you measure blood pressure regularly? | 53 |
| E | Technical information | |
| | 1. Error messages | 55 |
| | 2. Customer service | 57 |
| | 3. Technical data | 57 |
| | 4. Original spare parts and accessories | 58 |
| | 5. Metrological inspection | 58 |
| | 6. Explanation of symbols | 59 |
| | 7. Disposal | 60 |

Table of contents

| | | |
|----------|-------------------------------|----|
| F | Maintaining the device | 61 |
| G | Warranty | 62 |

Thank you for choosing the visomat® handy express wrist blood pressure monitor (hereinafter also referred to as the device).

These instructions for use are intended to help the user operate the device safely and efficiently. The unit must be used in accordance with the procedures contained in these instructions for use and must not be used for any other purposes.

It is important that you read all the instructions carefully before using the device. Please pay particular attention to the chapter "Important usage information" on page 42.

A**Intended use**

This device is intended for persons aged 12 or over for non-invasive self-measurement of the systolic and diastolic blood pressure, determining the pulse rate and calculating the pulse pressure.

This device uses the oscillometric method for measuring blood pressure and pulse rate.

Safety instructions**B****1. Important patient information**

- The device may only be used at the wrist and not on any other extremities.
- Taking blood pressure measurements on children requires specialist knowledge. Consult your doctor if you want to measure a child's blood pressure. Under no circumstances must the device be used on an infant.
- During pregnancy and in cases of pre-eclampsia, the device may only be used after consulting your doctor.
- Do not under any circumstances place the cuff over a critical area, e.g. a wound, aneurysm, etc. or on an arm with an arteriovenous shunt. Risk of injury! Any supply via an intravascular access point (infusion) or other medical monitoring devices could possibly be interrupted.
- Do not use the device without first consulting your doctor if you are having dialysis treatment or are taking anticoagulants, platelet aggregation inhibitors or steroids. These circumstances may cause internal bleeding.
- The device is not an instrument for the diagnosis of atrial flutter or atrial fibrillation. It shows only the possibility of the presence of such a condition. If the display indicating a cardiac arrhythmia appears (see page 49), it is essential to discuss this with your doctor.
- The display of the pulse frequency is not suitable for checking the frequency of cardiac pacemakers. Cardiac pacemakers and blood pressure monitors do not influence each other in terms of their mode of operation.
- The device is not suitable for use with electrosurgical devices.

- Do not use the device in explosive environments such as in the vicinity of flammable anaesthetics or in an oxygen chamber.
- Please wait a few minutes between successive measurements, otherwise the blood flow in the arm is interrupted for too long and injuries could occur.
- The device contains small parts that could be swallowed by children. Do not leave the device unattended with children.
- Do not take a blood pressure measurement on the same side as a mastectomy with simultaneous removal of the axilla lymph nodes.
- Discuss the measured values with your doctor. Do not assess the measurement results yourself. **Never change the dose of medicines prescribed by your doctor.**
- Before carrying out your self-measurement, please pay attention to the chapter “Important usage information” on page 42.

2. Important technical information

- The unit contains sensitive parts and must be protected from extreme temperature fluctuations, humidity, shocks, dust and direct sunlight.
- Do not knock or drop the device.
- The device is not waterproof.
- Do not apply pressure to the display. Do not place the device with the display facing downwards.
- The unit must only be operated with the cuff and accessories provided (see “Original spare parts and accessories” on page

58). Using different cuffs or accessories can lead to incorrect measurement results. The warranty will be invalidated if the unit is damaged by third-party accessories!


- A consistently good power supply to the device is necessary for fault-free blood pressure measurements.
 - Only use long-life alkaline batteries (LR03).
 - Always replace all batteries simultaneously during battery replacement. Do not mix new and old batteries or batteries of different types.
 - You will require 2 x 1.5 V batteries. Rechargeable batteries only have a voltage of 1.2 V and are therefore not suitable.
- If the device is not to be operated for an extended period of time, please remove the batteries. In principle, all batteries can leak.
- Do not remove batteries from the device when the device is switched on.
- The use of this device in the vicinity of mobile phones, microwaves or other devices with strong electromagnetic fields can cause malfunctions.
- Never open or modify the device or the cuff (except for replacing the batteries). If the unit has been opened, it must be subjected to a metrological inspection by an authorised institution.
- To avoid inaccurate readings, please adhere to the intended operating and storage conditions. See “Technical data” on page 57.
- The inflating and measuring procedure can be interrupted by pressing the Start/Stop button or by removing the cuff. The device then stops the inflation procedure and deflates the cuff.

1. Device description

- 1 Display
- 2 Memory button user 1 (memory recall button)
- 3 Memory button user 2 (memory recall button)
- 4 Start/Stop button
- 5 Cuff
- 6 Battery compartment

2. Display

- 1 Battery indicator
- 2 Motion control
- 3 Pulse signal display or irregular pulses
- 4 SYS = systole
- 5 DIA = diastole
- 6 PUL 1/min = pulse
Calculated pulse frequency per minute
- 7 User
- 8 Memory space or date/time
- 9 Pulse pressure control
- 10 Cuff position check
- 11 WHO classification (page 48)

3. Important usage information

The measurement results of automated blood pressure monitors can be influenced by the measuring location, posture, previous exertion and general physical condition. Please pay attention to the usage information in order to obtain accurate readings.

- Do not consume alcohol, caffeine or smoke for at least one hour before measuring.
- Rest for at least 5 minutes before taking the measurement. Depending on the previous degree of exertion, you may even need to rest for up to an hour.
- Expose the wrist. Under no circumstances may clothes hinder the blood flow to or from the arm, as this affects the blood pressure at the measuring point and can lead to incorrect readings.
- Do not move or talk during measuring.
- Breathe slowly and deeply. Do not hold your breath.
- Check the display for irregular pulse signals. If required, repeat the measurement when conditions have improved.
- In order to exclude side differences and to obtain comparable measurement results, it is important to always take the measurement on the same arm. Ask your doctor which side it is advisable for you to measure.
- Always take your daily blood pressure measurements at around the same time.
- Blood pressure is not a fixed value. It can go up or down by more than 20 mmHg in patients within a few minutes.

4. Initial operation of the device

Insert the batteries provided into the device.

The device has a date/time function that can be activated as required (see page 43). This function is deactivated on delivery.

When the date/time function is activated, the measurement results are saved along with their date and time. This can make it easier to view the readings at a later point.

5. Inserting/replacing batteries, battery icon

- Remove the lid of the battery compartment on the rear of the device by folding it up towards the cuff.
- Remove the old batteries from the device and insert the new ones. Please ensure correct polarity (marking in battery compartment).
- Close the battery compartment by clipping the battery cover back into the device.

When the “empty battery ” symbol first appears in the display, you are still able to take around 30 measurements. Please replace the batteries within this time frame.

6. Activating and setting the date/time

Activating and deactivating the date/time function

Press the Start/Stop button to activate the date/time function. Then immediately press one of the two Memory buttons while the full display is still visible. The function can be deactivated at any time using the same procedure.

Operating the device

When the date/time function is activated, the device continues to display the time even when it is switched off. The power consumption for this is extremely low.

Setting the date/time

Initially, the display will show the four-digit year number (Figure 1). Press Memory button 2 to reduce the year number or Memory button 1 to increase the year number. Confirm the entry by pressing the Start/Stop button.

Figure 1

You are now able to set the month. (Figure 2). Press Memory button 2 to reduce the value or Memory button 1 to increase the value. Set the correct month and confirm the entry by pressing the Start/Stop button.

Figure 2

Repeat these steps to enter the day (Figure 2) and the hours and the minutes (Figure 3).

Figure 3

After replacing the batteries, you will have to reset the date and time.

Changing the date/time


With the device switched off, hold down both memory buttons simultaneously for three seconds to enable the date and time to be set. Perform the setting as described on page 43.

7. Attaching the cuff

The special double-preformed comfort cuff is easy to attach and ensures a perfect fit during measurement. The materials used are tested for skin compatibility according to ISO 10993.

Operating the device

- Clear the wrist from all clothing or watches/jewellery since this may disturb the blood flow in the arm as it may impair the blood pressure at the place of measurement.
- Slide the device onto the centre of the inside of the forearm, directly on the wrist. The palm of your hand should be facing upwards. The distance between the wrist and the upper edge of the cuff should be 5–10 mm.
- Fasten the cuff using the hook and loop fastening.


The shape of the cuff is optimized to sit very tightly in order to avoid possible measurement errors.

8. Body position during the measurement

The measurements should be performed while seated. In special cases, measurements in a lying position are also possible.


Measurements while seated

Sit at a dining table, lean your back against the backrest of the chair and bend your arm so that your arm is exactly at heart level. Place your feet on the ground and do not cross your legs.

Measurements while lying down

Lay on your back, place the device on your wrist and at the height of your heart, e.g. by supporting it with a pillow. To relax before the measurement, take 2-3 deep breaths.

Please pay attention to the following:

- Do not bend your wrist.
- Keep your hand relaxed, not tightened, do not make a fist.
- It is important to keep your arm calm and not to speak during the measurement.

If the cuff is in a higher or lower position than the heart, the measured value might be higher or lower, respectively.

The results for measurements taken in a sitting and lying position may differ slightly.

9. Selecting the user

The device is designed for two users. The values for the two users are stored separately and can also be evaluated separately.

You may also use the two user memories to evaluate measurements taken in morning and those taken in the evening separately. To do this, carry out the measurement in the morning using user 1 and in the evening using user 2.

- Before taking the measurement, check that the correct user is displayed.
- To change the user, press the user button for whose memory bank you want to switch over to.
- Now carry out the measurement.

Also, you will be given another opportunity to assign the measurement to the correct user after the measurement has been taken (see page 48).

10. Measuring blood pressure

- Switch the device on using the Start/Stop button. The following display appears (Figure 1):


Figure 1

- After the device has completed its calibration against the ambient air pressure (Figure 2), the automatically controlled inflation process starts after three short beeps.


Figure 2

- Shortly after the inflation procedure begins, a check is carried out to ensure that the cuff has been applied tightly enough. A cuff that has been applied too loosely can result in inaccurate readings.

If the cuff is tight enough, the  symbol appears in the display (see Figure 3).

If the symbol is not displayed by the time the pulse signal starts, the cuff is too loose. In this case, you should repeat the measurement after you have applied the cuff correctly (see page 44).

- The measurement is carried out during inflation. The  symbol flashes in the display (Figure 3). In addition, beeps denote the pulse.


Figure 3

- Inflation stops as soon as the measured values are determined; the cuff is then automatically deflated.

- A long beep denotes the end of measurement.

- The determined systole, diastole and pulse values are shown in the display (Figure 4).

Figure 4

Operating the device

- When the measurement is finished, you will be given another opportunity to change the memory to which the measured value is to be saved. To do this, press the Memory button of the user for whom the measured value is to be saved.
- The device switches itself off after approximately 3 minutes. You can also switch the device off before this by pressing the Start/Stop button.

11. Classification of the measured values (WHO)

The device categorises the blood pressure readings according to the 1999 guidelines of the World Health Organisation (WHO). You can determine the classification after each measurement using the triangle in the display and the colour fields next to the display.

| Classification | Systolic pressure = top value mmHg | Diastolic pressure = bottom value mmHg |
|----------------------|--|--|
| Stage 3 hypertension | ≥ 180 | ≥ 110 |
| Stage 2 hypertension | 160-179 | 100-109 |
| Stage 1 hypertension | 140-159 | 90-99 |
| High-normal | 130-139 | 85-89 |
| Normal | 120-129 | 80-84 |
| Optimal | < 120 | < 80 |

The WHO classification is stored along with the measured values and both can be recalled from the memory.

Operating the device

Depending on age, weight and general health, the blood pressure values can be different. Only a doctor can determine the right blood pressure range for you and assess whether your blood pressure has reached a dangerous level for you. Discuss your blood pressure values with your doctor. **Never change the dose of medicines prescribed by your doctor!**

12. Pulse pressure control

The pulse pressure – not to be confused with the pulse beat – provides an indication of the elasticity of the blood vessels. A stiff system of vessels can have a negative effect on the cardiovascular system.

The heart works in two phases, the contraction phase (expulsion phase, systole) and the relaxation phase (filling or resting phase, diastole). The pressure difference between systole and diastole is called the pulse pressure or pulse amplitude. The higher the pulse pressure, the stiffer your vessels are thought to be.

Studies show that the cardiovascular risk increases if the pulse pressure value is permanently higher than 65 mmHg.

The device calculates the pulse pressure at the end of the measurement. If the value is greater than 65 mmHg, this is indicated by the  symbol on the display.

If your pulse pressure is permanently above 65 mmHg, consult your doctor.

13. Cardiac arrhythmia/atrial fibrillation

If the device detects irregular pulses during the measurement, the  symbol is displayed once the measurement is complete.

This can be triggered by cardiac arrhythmias, movement, talking or even by deep breathing.

The symbol is stored with the respective measurement.

 Symbol is flashing: The device has detected strong pulse disturbances with a very fast pulse. This may be **atrial flutter or atrial fibrillation**.

 Symbol is not flashing: Cardiac arrhythmia or irregular pulse detected.

The device is not an instrument for the diagnosis of atrial flutter or atrial fibrillation. It shows only the possibility of the presence of such a condition. **If this display should appear, it is essential to discuss this with your doctor to receive an accurate diagnosis.**

Measurement results accompanied by irregular pulses should be viewed as critical and repeated under more favourable conditions.

Caution: Do not confuse the irregular pulse display  with the flashing pulse display  that indicates the heart rate during every measurement (see page 51).

14. Motion control

Movements during the measurement may result in incorrect measurement values. The motion control registers stronger movements and these are indicated by the  symbol in the display. Other disruptive factors such as speaking, coughing or cuff noises can cause the symbol to be displayed. If this symbol appears, please repeat the measurement and avoid the disruptive factors.

The symbol is stored with the respective measurement.

15. Using the memory

The device stores the measured values in two memories (one per user). Each memory can store up to 120 readings and the average value (A). The average value is calculated from the last three measurements taken.

The measured results are automatically saved in the memory. If more than 120 measured values are saved, the oldest measured value (no. 120) is deleted so that the latest value (no. 1) can be recorded.

Data retrieval

Set the user whose data you wish to retrieve.

To retrieve data, press the Memory button of the user you wish to retrieve the data for while the device is switched off. The average value of the last three measurements is displayed along with an "A" (Figure 1).


Figure 1

Press the Memory button again to display the most recently measured values (Figure 2). If the date/time function is activated, the memory number (Figure 2), the date (Figure 3) and the time (Figure 4) are displayed alternately.


Figure 2

Figure 3

Figure 4

Press the Memory button repeatedly to recall further measured values.


You can switch between the two memories by pressing Memory button 1 or 2.

The stored data remain visible in the display for around 30 seconds. The device then switches itself off.


Deleting data

To delete individual measured values, recall the measured value you wish to delete.

Press the Memory button again and hold down until the display firstly flashes and is then replaced with dashes.

To delete the entire memory, retrieve the respective average value (A), press the Memory button again and hold down until the average value firstly flashes and is then replaced with dashes.

1. The systolic and diastolic blood pressure value

The circulatory system has the important task of supplying all organs and tissues in the body with sufficient blood, and of removing metabolites. The heart does this by contracting and expanding at a regular rhythm of around 60-80 times per minute. The pressure created by the blood flow on the artery walls when your heart contracts is referred to as systole. The pressure in the subsequent relaxation phase when the heart fills with blood again is referred to as diastole. In your daily measurement, you will measure both values.

2. Why you measure different values

Our blood pressure responds to internal and external influences like a sensitive measuring instrument. Affected by mental, physical and environmental influences, it varies all the time and never remains constant.

Reasons for fluctuating blood pressure values can be: moving, speaking, eating, consuming alcohol or nicotine, nervousness, inner tension, emotions, room temperature, recent urination or bowel movements, environmental influences such as movements and sounds, etc. Even changes in the weather and climate can have an effect on your blood pressure.

This explains why values measured at the doctor are often higher than those you obtain at home in your usual environment.

3. Why should you measure blood pressure regularly?

Even the time of day has an effect on your blood pressure. The values are at their highest during the day. In the course of the afternoon and in the evening, they drop slightly. They are low while you are sleeping, but rise again relatively quickly once you get up.

D What you should know about blood pressure

One-off and irregular measurements therefore say little about your actual blood pressure. A reliable assessment is only possible when measurements are taken regularly. Discuss the measured values with your doctor.

Technical information

E

1. Error messages

| Error | Possible cause | Remedy |
|---|--|--|
| Display Err - 1 | Measurement error; measurement could not be carried out due to movement or speaking. | - Repeat the measurement - Do not move your arm - Do not speak |
| Display Err - 2 | Inflation error; pressure built up too quickly or too slowly. | - Repeat the measurement - Do not move your arm - Do not speak - Check wrist circumference. Call the Service hotline if this occurs repeatedly. |
| Display Err - 3 | Excessive cuff pressure. The cuff may have been inflated to the maximum by moving the arm or the body. | - Repeat the measurement - Do not move your arm - Do not speak |
| Display 0 0 | Start/Stop button was accidentally pressed while the batteries were being replaced. | Switch unit off and then on again using the Start/Stop button. |
|  | Irregular pulse, possible cardiac arrhythmia. | Repeat measurement after resting for 3-5 minutes. Listen out for the device's pulse beep signals. Possible influences on the measurement results (depending on the severity of the arrhythmia) must be discussed with your doctor! |
|  | Possible atrial flutter or atrial fibrillation. | Repeat measurement after resting for 3-5 minutes. Listen out for the device's pulse beep signals. Possible influences on the measurement results (depending on the severity of the arrhythmia) must be discussed with your doctor! |
|  | The motion control detected movement during the measurement. | - Repeat the measurement - Do not move your arm - Do not speak |

| Error | Possible cause | Remedy |
|--|--|---|
| Measured values/pulse values are extremely high or low. | Too little rest before taking the measurement. Movement or talking during the measurement. | Repeat the measurement after a break of 3-5 minutes. Pay attention to the usage information on page 42. |
| | Incorrect measurement posture | Repeat the measurement and maintain the correct posture. (Page 45) |
| Measurements taken in direct succession produce different measured values. | Blood pressure is not a fixed value. It can go up or down by more than 20 mmHg within a few minutes. | In order to obtain comparable results, always carry out the measurements under the same conditions. |
| Measurement was interrupted, the cuff deflated to 0 and the measurement restarted. | Device detects a fault or that the diastolic value is too low. | Device shows no fault. Repeat measurement |
| | Movement during measurement | If necessary, interrupt measurement and repeat after 5 minutes' rest. |
| Dashes in display after pressing Memory button | No measured values have been stored. | Repeat measurement |
| The measurement is interrupted and the  symbol appears. | Batteries flat | Replace batteries. |
| The  symbol appears during measurement but the measurement continues. | Low batteries | You can still perform approx. 30 measurements, but should keep spare batteries available. |
| No display after device is switched on | Batteries inserted incorrectly | Check position of batteries. |
| | Batteries flat | Replace batteries. |
| | Battery contacts dirty. | Clean battery contacts with a dry cloth. |
| Date and time are not visible in the display | Function was unintentionally deactivated | Activate the date and time function as described on page 43. |

2. Customer service

Device repairs may only be performed by the manufacturer or an expressly authorised body. Please contact:

UEBE Medical GmbH
 Zum Ottersberg 9
 97877 Wertheim, Germany
 Tel.: +49 (0) 9342/924040
 Fax.: +49 (0) 9342/924080
 Email: info@uebe.com
 Website: www.uebe.com

3. Technical data

| | |
|--|--|
| Model: | visomat® handy express |
| Size: | 88 mm x 64 mm x 26.4 mm |
| Weight: | approx. 101 g without batteries and cuff |
| Display: | LCD display (liquid crystal display) for measured values and check displays |
| Memory: | 2 memories each capable of storing 120 readings |
| Measurement method: | Oscillometric determination of systole, diastole and pulse |
| Reference method of the clinical test: | Auscultatory measurement |
| Pressure display range: | 0-300 mmHg |
| Measurement range: | Systolic: 50-250 mmHg Diastolic: 40-150 mmHg Pulse measurement: 30-199 pulses/minute |
| Measurement precision: | Blood pressure measurement (cuff pressure): ± 3 mmHg Pulse rate: ± 5% |
| Power supply: | 2 x 1.5 V AAA LR03 alkaline manganese batteries |
| Power consumption: | 2 W |

| | |
|---------------------------------------|---|
| Cuff: | Cuff type H7, wrist circumference 12.5 cm - 21.5 cm |
| Operating conditions: | Ambient temperature 10°C to 40°C, relative humidity 15% to 85% non-condensing, air pressure 700 to 1060 hPa |
| Conditions for storage and transport: | Ambient temperature -20°C to +60°C Relative humidity 15% to 95%, non-condensing |
| Inflation and deflation: | Regulated automatically |
| IP rating: | IP 20: Protected against solid foreign bodies with a diameter from 12.5 mm, no protection against water. |
| Protection against electric shock: | Internal power supply, applied part type BF (cuff) |
| Operating mode: | Continuous operation |
| Classification: | Internal power supply using batteries. |

4. Original spare parts and accessories

The following original spare parts or accessories are available through specialist retailers:

- Cuff type H7 (for wrist circumference 12.5-21.5 cm)
Part no. 2105601, PZN-13975619

Subject to technical modifications.

5. Metrological inspection

Generally, it is recommended that a metrological inspection is performed every two years. However, professional users in Germany are obligated to do so in accordance with the "Medical Products Operator Ordinance" (Medizinprodukte-Betreiberverordnung). This can be performed either by UEBE Medical GmbH, an authority responsible for metrology or an authorised maintenance service. Please refer to your national regulations.

Upon request, responsible authorities or authorised maintenance services receive a "Test instruction for metrological inspection" from the manufacturer.

Please only submit the device for metrological inspection together with the cuff and the instructions for use.

Important: No modifications, e.g. opening the device (except to replace the batteries), may be made to this device without the manufacturer's permission.

6. Explanation of symbols

 This product complies with Council Directive 93/42/EEC of 5 September 2007 concerning medical devices and bears the mark CE 0123 (TÜV SÜD Product Service GmbH).

Degree of protection against electric shock TYPE BF

Consult the instructions for use

Keep dry

Device serial number

Reference number = item number

IP20 Protected against solid foreign bodies with a diameter from 12.5 mm, no protection against water.

Manufacturer

7. Disposal

Technical equipment and batteries do not belong in household waste. They must be disposed of at appropriate collection and disposal points.

- Keep the device clean. Check for cleanliness after use. Please use a soft, dry cloth for cleaning. Do not use benzene, thinners or other strong solvents.
- The cuff may absorb sweat and other fluids. Inspect the cuff for marks or discolourations after each use. Please use a synthetic detergent for cleaning and carefully wipe the surface. Do not scrub or machine-wash. Air-dry thoroughly.
- Make sure when storing the unit that no heavy objects are left lying on the unit.
- Do not apply force when folding or bending the cuff.
- If the device was stored at temperatures below 0 °C, allow it to stand at room temperature for at least 1 hour before using it.
- If the device is not to be operated for an extended period of time, please remove the batteries. Store the batteries out of the reach of children.

The device has been manufactured and tested with all due care. Nevertheless, in the event of defects upon delivery, we provide a warranty under the following conditions:

During the warranty period of 3 years from the date of purchase, we remedy such defects at our discretion and at our expense in our plant through repair work or replacement delivery of a defect-free device.

The warranty does not cover normal wear and tear of wearing parts or damage caused by failure to observe the instructions for use, improper handling (e.g. unsuitable power source, breakage, leaking batteries) and/or disassembly of the device by the buyer. Furthermore, the warranty does not constitute grounds for asserting claims for damages against us.

Warranty claims can only be asserted during the warranty period and upon presentation of the proof of purchase. In the case of a warranty claim, the device along with the proof of purchase and a description of the complaint must be sent to:

UEBE Medical GmbH
Service-Center
Zum Ottersberg 9
97877 Wertheim, Germany

The cost of returning the device to our factory shall be borne by the sender. Complaints that are sent back without prepayment will not be accepted by UEBE.

The statutory claims and rights of the buyer against the seller (for example, defect claims, producer liability) are not restricted by this warranty.

Please note: In the case of a warranty claim, please make sure to enclose the proof of purchase.

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© Copyright 2017 UEBE Medical GmbH