




 **BATEC**
MOBILITY



















BATEC • ELECTRIC



The BATEC ELECTRIC offers you a 2 in 1 approach to your needs: a unique solution that will revolutionise your mobility outdoors, without having to renounce to the advantages of your manual wheelchair. It is the definitive mobility tool.



SPECIFICATIONS

- MATERIAL: 7005 Aluminium
- MOTOR: Brushless, 500 W and 36 V
- BATTERY: LIFEPO4, 9 Ah Lithium, without memory effect (25-30 km) with wireless automatic connection system
- DOUBLE USB PORT for charging a smartphone or a tablet
- Ultra-quick CHARGER (4 hours)
- TYRE: 20" Maxxis ultra-gripping (54-406)
- WHEEL RIMS: 20" double wall aluminium
- BRAKES: Double Shimano 203 mm mechanical disc brake with Shimano Deore levers
- HEADSET: FSA Orbit MX
- HANDLEBAR: FSA V-Drive Lowriser
- BATEC EASY-FIX anchor system (Pat. Pend.)
- BATEC EASY-FIT regulation system (Pat. Pend.)
- BATEC SAFE-FEET stand frame (Pat. Pend.)
- FRONT LED LIGHT: 900 lumen, connected to the main battery
- DOUBLE REAR LED LIGHT on the stand frame
- CONTROL PANEL with battery indicator and power regulator
- CYCLE COMPUTER with LCD screen and 4 functions
- REMOVABLE QR COUNTERBALANCE WEIGHTS to improve traction
- WEIGHT: 14.9 Kg (+ 4 Kg battery + 5 Kg counterbalance weights)

MEASURES IN MM



COLORS



BATEC • MANUAL



Connect a BATEC MANUAL to your wheelchair and discover a new way of moving about and exercising.



SPECIFICATIONS

- MATERIAL: 7005 Aluminium
- TYRE: 20" Maxxis ultra-gripping (54-406)
- WHEEL RIMS: 20" double wall aluminium
- SHIFTER: Sram Dual Drive with 24 speeds (3 speed internal hub + 8 gears)
- BRAKES: Double brake with 180 mm Shimano mechanical disc brake with Shimano Deore levers
- HEADSET: FSA Orbit MX
- CRANKSET: Shimano 170 mm aluminium
- BATEC EASY-FIX anchor system (Pat. Pend.)
- BATEC EASY-FIT regulation system (Pat. Pend.)
- BATEC SAFE-FEET stand frame (Pat. Pend.)
- FRONT LIGHT with 7 LEDS with 4 functions
- REAR LIGHT with 2 LEDS with 2 functions
- CYCLE COMPUTER with LCD screen and 8 functions
- REMOVABLE QR COUNTERBALANCE WEIGHTS to improve traction
- WEIGHT: 13 Kg (+ 5 Kg counterbalance weights)

MEASURES IN MM



COLORS



BATEC • HIBRID



Thanks to its electric assisted motor, the BATEC HIBRID will allow to you enjoy healthy physical exercise.



SPECIFICATIONS

- MATERIAL: 7005 Aluminium
- MOTOR: Brushless, 500 W and 36 V
- BATTERY: LIFEPO4, 9 Ah Lithium, without memory effect (25-30 km) with wireless automatic connection system
- DOUBLE USB PORT for charging a smartphone or a tablet
- Ultra-quick CHARGER (4 hours)
- TYRE: 20" Maxxis ultra-gripping (54-406)
- WHEEL RIMS: 20" double wall aluminium
- SHIFTER: Shimano Alivio 8 speed
- BRAKES: Double Shimano 203 mm mechanical disc brake with Shimano Deore levers
- CRANKSET: Shimano 170 mm aluminium
- BATEC EASY-FIX anchor system (Pat. Pend.)
- BATEC EASY-FIT regulation system (Pat. Pend.)
- BATEC SAFE-FEET stand frame (Pat. Pend.)
- FRONT LED LIGHT: 900 lumen, connected to the main battery
- DOUBLE REAR LED LIGHT on the stand frame
- CONTROL PANEL with battery indicator and power regulator
- CYCLE COMPUTER with LCD screen and 8 functions
- REMOVABLE QR COUNTERBALANCE WEIGHTS to improve traction
- WEIGHT: 19,9 Kg (+ 4 Kg battery + 5 Kg counterbalance weights)

MEASURES IN MM



COLORS



Social project

Batec Mobility is a socially-oriented company created by people with disabilities with a view to making a mark on this very sector of society, by making a positive social impact.

Our greatest impact is on **Batec** users, whose life we manage to change by improving their mobility. We are proud that more than 94 % of our users state that **Batec** has substantially changed their lives for the better.

Secondly, **Batec Mobility** also creates jobs for disabled people. Lastly, we collaborate closely with institutions, promoting an active lifestyle for people with disabilities.





Batec Mobility is based in Barcelona. The nature of this Mediterranean city, which invites people to enjoy the open air in good company, inspired us to create Batec, to enable people with physical disabilities from across the world to lead an active, dynamic lifestyle. Barcelona is also a large city with an old tradition in industrial design, which is recognised worldwide for its added value and quality. This was also decisive when it came to developing our products, which are essentially useful but also very attractive to the eye.



© Batec Mobility S. L., 2013
www.batec-mobility.com

Models: Alba, Pablo, Javier, Patricia and Clara
Batec users: Ana, Axel, Alina, Joan, Chus, Teo, Ana, Oriol, Sofia, Sergio, José Miguel, Bàrbara, Rafael, Manuel, Hugo, Rafael, Òscar and Maria Andrea
Staff: Òscar, Pau, Javier and Carlos
Institution: Aspaym Catalunya

Pictures of the catalogue by Asensi Carricondo
www.asenicarricondo.net
Design and layout by Cristina Serrano
www.cristinaserrano.net

Translation by Duthie Translations

Printed in Spain by Mediaactive S. L.



Created by and for disabled people, **Batec Mobility** revolutionises your mobility

www.batec-mobility.com

